



People.Health.Care.

VENOSAN® 6000

Classic compression stockings.



VENOSAN® 6000


Classic compression stockings.

The economical stocking which meets all medical-grade graduated compression requirements.




Classic appearance

Slight transparent construction for everyday wear.



Simple construction

Designed with common highly elastic yarns for an affordable option that meets clinical needs.



Anatomically-knit

Provides maximum fit and comfort which stay supported on the leg.



Easy care

Machine wash and air dry naturally overnight.



Swiss comfort

Manufactured in Switzerland to the highest standards.



Available Styles



VENOSAN® 6000 medical compression stockings offer mid-containment of Stage 1–2 chronic venous insufficiency, varicose veins, varicoses during pregnancy and support post-vein surgery. Its classic yarns and simple construction make it an affordable everyday solution, while also ensuring it meets the therapeutic compression needs patients deserve.

Indications

- Chronic Venous Insufficiency (CVI) (Stage 1–2)
- Varicose veins
- Varicosis during pregnancy
- Post vein surgery

Options

Compression 20–30 mmHg, 30–40 mmHg

Toe Style Open toe, Closed toe

Fit Standard

Colours

Beige

Black



People.Health.Care.


VENOSAN® 6000
Classic compression stockings.

**Circumference measurements |
Mesures de circonférence (cm)**

	SMALL	MEDIUM	LARGE	X-LARGE
cB	20-22	23-25	26-28	29-31
cC	30-39	33-41	36-44	39-47
cD	29-38	32-40	35-43	38-46
cG	49-63	54-66	58-70	63-74

Length | Longueur (cm)

		short	long
AD	ℓD	34-39	39-44
AGH	ℓG	59-67	68-78
AT	ℓK	67-74	74-85



Material content | Composition
70% polyamide, 30% elastane



L&R CAN INC. · 2820 Argentia Road, Unit 2, Mississauga, ON L5N 8G4 · 888-836-6726 · Fax: 888-434-5347 · inquiries@ca.LRmed.com · www.Lohmann-Rauscher.ca

VENOSAN® is a registered trademark of SWISSLASTIC AG, ST. GALLEN, Switzerland — an L&R Company