




People.Health.Care.

VENOSAN[®] 4000

High-performance support.



www.Lohmann-Rauscher.ca

VENOSAN® 4000


Climate effect compression stockings and armsleeves.

The construction of the microfiber mesh in the VENOSAN® 4000 compression stocking ensures lightness and well-being.



Lightness and breathability

Air-permeable yarns and knitting construction offer maximum comfort by wicking away moisture.



Microfiber construction

The velvety touch of the high tech TACTEL® microfibres — many times thinner than a human hair — offers outstanding comfort.



Fashionably elegant

Wearable in all seasons. It has a contemporary stylish look, thanks to the semi-sheer transparency.



Supports self-care

Easy to put on and take off with its stretchable construction.



Swiss comfort

Manufactured in Switzerland to the highest standards.




Available Styles



VENOSAN® 4000 medical grade compression stockings and arm sleeves offer a climate effect microfiber mesh construction for mild to moderate venous edema and preventative support. It is a reliable, everyday companion with excellent wearer satisfaction, because of its velvety touch and wickability with microfiber yarns providing maximum comfort and thereby a better quality of life.

Indications

- Chronic Venous Insufficiency (CVI) (Stage 1)
- Pregnancy-induced varicose veins
- Tired, heavy and aching legs
- Post vein surgery



Options

Compression 20–30 mmHg, 30–40 mmHg

Toe Style Open toe, Closed toe

Fit Standard, Made to Measure

Colours



A TACTEL® microfibre yarn

B Human hair



People.Health.Care.

VENOSAN® 4000


Climate effect compression stockings and armsleeves.

Measurements | Mesures de circonférence (cm)

	SMALL	MEDIUM	LARGE	X-LARGE
cB	20-23	23-26	26-29	29-31
cC	29-39	33-42	36-45	39-47
cD	28-38	32-41	35-44	38-46
cG	48-63	54-67	58-72	63-74

Length | Longueur (cm)

		short	long
AD/ADH	ℓD	35-40	37-42
AGH	ℓG	63-72	68-78
AGG/AT/ATU	ℓK	67-74	74-85




Measurements | Mesures de circonférence (cm)

	SMALL	MEDIUM	LARGE	X-LARGE
cC	14-17	17-20	20-23	23-26
cD	22-26	26-30	30-34	34-38
cF	24-30	30-36	36-42	42-48

Length | Longueur (cm)

		long
CG	ℓCG	40-50



Material content
71% TACTEL® polyamide, 29% elastane



L&R CAN INC. · 2820 Argenta Road, Unit 2, Mississauga, ON L5N 8G4 · 888-836-6726 · Fax: 888-434-5347 · inquiries@ca.LRmed.com · www.Lohmann-Rauscher.ca

VENOSAN® is a registered trademark of SWISSLASTIC AG, ST. GALLEN, Switzerland — an L&R Company
TACTEL® fibre is a registered trademark of INVISTA and is used under license to NILIT for nylon fibre products.